


KANGJI





LIGATION SOLUTION TITANIUM CLIPS

Easy Loading with Precision

Titanium Ligation Clips

- Made of titanium with easy loading system
- With inner serrations to ensure secure ligation
- Compatibility with MRI




Right Applicator Jaw

Keep jaw angle 90° to the cartridge

Product	Size	Size Range of Vessel / Tissue Bundle	Code	Clips/Cartridge	Cartridge/Box
Titanium Ligation Clips	S (LT100) 	0.2-2.9mm	111Y.204	6	20
	M (LT200) 	0.2-4.0mm	111Y.203	6	20
	ML (LT300) 	0.4-5.5mm	111Y.202	6	20
	L (LT400) 	0.9-6.5mm	111Y.201	6	20

Titanium Clip Applier single Action




Product	Code	Clip Size	Type	Shaft Diameter (MM)	Working Length (MM)	Optional Working Length(MM)
Titanium Clip Applier	104.001B	S(LT100)	Endo-surgery, single action	Φ5	300	360
	104.001A	M(LT200)	Endo-surgery, single action	Φ5	300	330
	101.001A	ML(LT300)	Endo-surgery, single action	Φ10	330	360
	101.001AF	ML(LT300)	Endo-surgery, single action, detachable	Φ10	330	/
	101.001B	L(LT400)	Endo-surgery, single action	Φ10	330	360
	101.001BF	L(LT400)	Endo-surgery, single action, detachable	Φ10	330	/

Titanium Clip Applier Dual Action




Product	Code	Clip Size	Type	Shaft Diameter (MM)	Working Length (MM)
Titanium Clip Applier	101.111A	ML (LT300)	Endo-surgery, dual action	Φ10	330
	101.111AF	ML (LT300)	Endo-surgery, dual action, detachable	Φ10	330
	101.111B	L (LT400)	Endo-surgery, dual action	Φ10	330
	101.111BF	L (LT400)	Endo-surgery, dual action, detachable	Φ10	330
	101.111ABF	ML/L	Endo-surgery, dual action, detachable	Φ10	330

Titanium Clip Applier 90°



Product	Code	Clip Size	Type	Shaft Diameter (MM)	Working Length (MM)
Titanium Clip Applier	101.002	ML (LT300)	Endo-surgery, 90°	Φ10	330/380
	101.002F	ML (LT300)	Endo-surgery, 90° detachable	Φ10	330

Titanium Clip Applier Open Surgery



Product	Code	Clip Size	Type	Working Length (MM)	Optional Working Length (MM)
Titanium Clip Applier	101.138CP	S (LT100)	Open Surgery, 15°	267	/
	101.138BP	M (LT200)	Open Surgery, 15°	267	230
	101.138AP	ML (LT300)	Open Surgery, 15°	267	/
	101.138P	L (LT400)	Open Surgery, 15°	267	200

KANGJI

Hangzhou Kangji Medical Instrument Co., Ltd.

Stock Code: 9997.HK

Address | No.1668 Chunjiang East Rd. Economic Development Zone,
Tonglu, Hangzhou, Zhejiang 311501, China.

Tel. | +86-571-69900029

Email | kangji@kangji.com

Website | <https://www.kangjimed.com>



Follow us for more